



**NEW!**



**custo watch**

**The Holter Monitor.**

## custo watch:

### The destination of a long journey.

**Our vision:** to make Holter ECG recordings as comfortable as possible for the patient.

**The interim stages:** a chest-attached device instead of a hip-mounted recorder, 3 electrodes instead of 7, very short patient cables instead of long ones.

**The result of further consistent development work:** custo watch. The 24-h ECG is recorded with a lightweight chest belt, transmitted wirelessly to the ECG watch and it is saved there. It comes without adhesive electrodes or wires but with a high degree of patient comfort and an attractive design.

custo watch provides noticeable and measurable benefits for both the patient and the doctor. For our part, we have our next goals in mind.



ECG marker button

ECG check before starting to record



## custo watch:

### The advantages of a sophisticated system.



custo belt 3

#### custo belt 3

- Holter ECG without adhesive electrodes
- Novel belt technology for 3-channel ECG recording
- Reusable
- Machine washable
- High wearing comfort
- Certificate of biocompatibility



custo guard 3

#### custo guard 3

- 3-channel ECG transmitter
- Waterproof, extremely lightweight
- Energy-saving wireless transmission of ECG signal to custo watch



custo watch

#### custo watch

- Signal testing of ECG channels
- Adjustable patient display
- Interchangeable bracelets
- Waterproof, extremely lightweight
- Recording time: 24 h



custo docking station

#### custo docking station


- Charging of custo guard and custo watch
- Fast data transfer of Holter ECG data to the PC
- Compatible with Holter ECG evaluation software custo tera

## custo diagnostic:

# One platform for all applications.

The consistent user interface of custo diagnostic manages all patients in proper cardiorespiratory diagnostic modules.

All required medical data are exchanged with the superordinate administrative system (HIS, EPR) via standardised interfaces (HL7, GDT, DICOM).



# Fax: +49 89 7109810

custo med GmbH • Leibnizstraße 7 • D-85521 Ottobrunn • Phone: +49 89 710 98 00 • Email: [info@customed.de](mailto:info@customed.de)

## Please send me a non-binding offer for custo watch

Kindly fill in the form and send it to custo med by telefax or by mail.

Please send me information on the subjects checked below:

☐ Resting/Stress Test ECG

☐ Pulmonary Function

☐ Holter ECG

☐ ABPM

☐ Cardiorespiratory Polysomnography

☐ Prevention and Rehabilitation

Address:

Name \_\_\_\_\_

Street \_\_\_\_\_

Postcode/Town \_\_\_\_\_

Phone \_\_\_\_\_

Telefax \_\_\_\_\_

Email \_\_\_\_\_

Doctor's office stamp